**EN**

POLITRAY is an acrylate for the fabrication of individual impression trays by using the cold curing.

**Indications:** individual impression tray

**DIRECTIONS FOR USE**

- Preparation: In order to cover the undercuts on the model wax is to be applied. The gypsum model is to be isolated by an isolating fluid.
- Mixing and fabrication of the individual impression tray: The recommended mixing ratio is 30 g powder : 10 ml liquid (ratio of volumes 3,5:1). Place the liquid into a mixing jar and add the powder little by little. In order to achieve a homogeneous mixture the liquid and the powder are to be mixed with a spatula for about 30 sec. The consistency of the mixture after that is suitable for the preparation of the impression tray (The mixture doesn't stick anymore and separates easily from the sides of the mixing jar). Knead the mixture, place it on the model and form the impression tray. Trim off all the excesses. The mixture at room temperature of 22°C can be handled for a period of 4,5 minutes. The working time is subject to the surrounding room temperature.
- Polymerisation: At room temperature of 22°C the autopolymerisation of the individual impression tray mixture takes place in 30 minutes.
- Trimming: After the polymerisation the impression tray is removed from the model and trimmed according to procedures common in dental laboratories

**Storage** Store at temperature 5-25°C. Avoid direct exposure to sun. Store in dry place.

**Warning** Close container immediately after use. Do not use material after expiry date. Keep out of reach of children. Avoid direct contact between skin and uncured material as well as inhalation of the liquid. In case of known allergies or hypersensitivity to acrylic materials do not use the product. In case of severely immuno compromised patient, the doctor should decide and weigh the benefit/risk with regards to the individual patient. Due to lack of information regarding the safety of PMMA exposure for pregnant women, we advise to avoid using the device for pregnant women. The final use in this special group of patients should be weighed by the dentist according to risk/benefits ratio for each patient.

In case that any serious incident has occurred in relation to the device it should be reported to the manufacturer (Polident d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, where the user is established.

**Warning of liquid componente** Liquid is highly flammable. Irritating to eyes, respiratory system and skin. May cause sensitization by skin contact. Keep container in a well ventilated place. Keep away from sources of ignition. NO smoking. Do not empty into drains. Take precautionary measures against static discharges.

**Packaging**

Politray powder 350g, 700g  
Politray liquid 150ml, 500ml

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