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INSTRUCTIONS FOR USE

1. Storage and working temperature

Powder and Liquid should be stored at room temperature (23 ±1°C) before using.

Storage and working temperature of investment powder and liquid are an important factor in determining the setting time and expansion, and hence the fit and the surface roughness of the pressed ceramic objects.

If the special liquid is stored at temperature below 5°C, it will be subjected to be frozen and will not be suitable to be used.

Working temperature: Use at **23±1°C** room temperature. Higher working temperatures reduce working time, lower temperatures prolong working time.

2. Preparations before investing

Polipress investment material is a very fluid investment, which can be used without any wetting agents. If a wetting agent is used, be sure to totally dry the surface prior to investing.

Use adequate silicone ring systems for the ceramic press technique in size 100g and 200g.

Powder/liquid ratio: Investment material powder is mixed with Polipress special liquid which is diluted with distilled water.

Exact powder/liquid measurement is necessary to obtain stable results.

| Ring size | Powder | Liquid |
|-----------|--------|--------|
| Small | 100 g | 22 ml |
| Large | 200 g | 44 ml |

Expansion: The expansion of the investment can be controlled by the quantity of the distilled water mixed with the special liquid. Higher is the concentration of the special liquid, higher is the total expansion of the investment. Use only distilled water to dilute.

The concentration can be freely adapted based on the working experience of the technician.

Polipress liquid dilution chart:

| Technique | Conc. (%) Polipress Liquid | 100 g ring Liquid : distilled H2O | 200 g ring Liquid : distilled H2O |
|--|----------------------------|-----------------------------------|-----------------------------------|
| Inlays | 40% | 9 ml : 13 ml | 18 ml : 26 ml |
| Partial crowns, veneers, single crowns | 45% | 10 ml : 12 ml | 20 ml : 24 ml |
| 3-unit bridges | 40% | 9 ml : 13 ml | 18 ml : 26 ml |
| Press on metal and ZrO2 frameworks | 50% | 11 ml : 11 ml | 22 ml : 22 ml |

Note, influences such as an increased room temperature, increased material temperature, model type (wax, plastic, wax/plastic combination) etc. may affect the processing properties or final result.

3. Mixing

Pre-mix powder and liquid thoroughly by hand with a spatula. Be sure all powder is wetted out by the liquid to give a uniform mixture.

Mix for 60 seconds under vacuum (400 rpm).

Always use clean mixing bowl and check vacuum level. Insufficient vacuum leads to differences in fit and air-bubbles of the pressed objects.

Working time: At least **4 minutes** pouring time including the mixing time at room temperature (23°C).

Working time/pouring time depends on the temperature of powder & liquid and on room temperature. Higher temperatures shorten the working time.

4. Investing

Invest under gentle vibration (low frequency).

Polipress investment material is very fluid, so strong vibration is not necessary and not advisable. From the moment the ring is totally filled, stop vibration immediately and do not touch the investment until set.

Setting time: Leave to set for 20 minutes from start of mixing.

Best results are obtained by putting immediately into a preheated furnace after 20 min.

Scrape the top surface of the investment ring with a sharp knife/sandpaper.

5. Heating up procedure

| | Quick heating | Conventional step heating |
|--|---------------|---------------------------|
| Insertion temperature | 850°C | 250°C |
| Holding time at 250°C | | 60 mins |
| 250°C ? 570°C, heating rate | | 9°C / min |
| Holding time at 570°C | | 30 mins |
| 570°C ? End temperature, heating rate | | 9°C / min |
| End temperature | 850 - 900°C | 850 - 900°C |
| Holding time at end temperature (depending on size of cylinder) | 45 - 60 mins | 45 - 60 mins |

Preheating temperature: Always place the set investment ring in the burnout furnace at 850°C.

When the mould is placed into the furnace, preheated furnace temperature should not be higher than 850°C. Then heat up to higher temperature (900°C) if required.

Due to aggressive burning out, do not open the furnace during the first 15 minutes.

When several investment rings are put into the furnace at the same time, prolong the heating period of each step with 10 minutes.

Holding time, after the preheating temperature of 850°C:

100 g ring – min. 45 minutes
200 g ring – min. 60 minutes

6. Pressing the ceramic and press plungers

Follow carefully the instructions for use of the pressable ceramic.

Press sequence should start as soon as possible after removing the investment ring from the furnace.

Standard Alox plungers and one-way investment plungers can be used.

7. Cooling

Follow carefully the instructions for use of pressable ceramic.

Safety recommendations and hazard warnings

1. Investment materials contain quartz. Do not inhale dust! Risk of lung damage (Silicosis, lung cancer). Wear protection mask Type FFP 2 – EN 149:2001.

Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl. Rinse the empty investment material bag with water before disposal.

2. To avoid the formation of dust when removing the investment material from the object, place the cooled mould into water for a short time.

3. When blasting the object, always use a fine-dust filter extraction system.

Expiry date: 2 years from the manufacturing date

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