EN-K30 S Heat-curing acrylic resin for full and partial dentures by using the **packing pressing** technique

General description: Denture base polymer material K-30 S is intended for the fabrication of full and partial dentures by using the packing/pressing technique and classical heat curing. It is suitable for edentulous/partial edentulous patients who need to restore functions of missing teeth (mastication, phonetics) and also to provide an aesthetic appearance for the patient.

The material is supplied in the form of liquid and powder that shall always be used together by a dental health care professional (dental technician).

Material characteristic	Specification	Method according to
Flexural strength	> 65 MPa	EN ISO 20795-1
Flexural modulus	> 2000 MPa	EN ISO 20795-1
Sorption	< 32 µg/mm³	EN ISO 20795-1
Solubility	< 1.6 µg/mm³	EN ISO 20795-1
Residual monomer	< 2.2 %	EN ISO 20795-1

Variants: Available in various colours and quantities. Consult your distributor.
Indication: Edentulism, partial edentulism

Composition:
Powder: PMMA-based powder, benzoyl peroxide, pigments
Liquid: Methyl methacrylate, trimethacrylate, dimethacrylate

Contraindications:
Do not use it in the case of known allergy/hypersensitivity to acrylates or any other component. A safety data sheet is available on the web page.
In the case that any serious incident has occurred in relation to the device, it should be reported to the manufacturer (Polident d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, where the user is established.

DIRECTIONS FOR USE:
Do not use K-30 S for direct relining of the denture base. Visually inspect the liquid. If you notice any sediment or gel, do not use the liquid. Discard it in accordance with local legislation.
We recommend temperature conditioning of the liquid and powder components at room temperature before use.

- Preparation: Flasking and eliminating the wax from the flask should be done following common dental laboratory procedures. Add liquid detergent to the boiling water and continue rinsing with boiling water until the wax is thoroughly removed from the teeth. Coat the surface of the plaster model with an isolating flacon and ensure it is thoroughly dried. To ensure a tight bonding between acrylic teeth and denture base plate, apply the K-30 S liquid on the roughened acrylic teeth. Mechanical retention is recommended and is essential when bonding teeth made of ceramic, composites, or other materials.

- Mixing and packing/pressing: The recommended mixing ratio is 23 g powder to 10 g of liquid (volume ratio 3:1). Pour the liquid into a mixing jar and gradually add the powder. Mix the liquid and powder with a spatula for about 30 sec. Cover the mixing jar and set it aside for around 30 minutes at room temperature (approx. 22 °C). After that, the consistency of the mixture is suitable for the packing (the mixture should no longer be sticky and should separate clearly from the walls of the mixing jar). The working time of the mixture is of approximately 15 minutes. The dough time as well as the working time is subject to the temperature of the liquid and the powder before mixing and the room temperature. Knead the mixture with your fingers and place it into the flask. The use of a cold flask is recommended. Close the flask carefully and put it into the press. Gradually increase the pressure in the press up to 80 bar and leave it at that level for 10 minutes. Insert the flask into a clamp.

- Polymerisation: You can choose among the following polymerization types:
Type 1: Immerse the flask in its clamp in cold water, heat the water to 100°C and polymerize for 35-45 minutes.

Type 2: Immerse the flask in its clamp in cold water, heat the water in 30 minutes up to 70°C, polymerize for 30-45 minutes. Then heat the water in 30 minutes up to 100°C, and polymerize for 30-45 minutes.

- Cooling: Cool down the flask for 30 minutes in the open and after that in a cold water.

- Trimming: Once the flask is completely cool, open it and remove the plaster model. Trim and polish the prosthesis according to common dental laboratory procedures.

Repairs: The use of the cold curing material Polirepar S is recommended, the colours of which matches the acrylate for the fabrication of prothesis by using the heat curing (K-30 S).

NOTE! Dental restorations made of acrylys always contain a small proportion of residual monomer. The residual monomer content is below 2,2% if the polymerisation is performed according to the instructions for use. For further free monomer reduction it is recommended to store a restoration in water for at least 24 hours. If this step is not followed, the patient to place the appliance in water for a period whenever it is not being worn.
Not following the instructions for use can lead to elevated residual monomer outside the permitted limits, poor physical properties or insufficient surface properties leading to biofilm formation.

Information to be provided to the dentist:
In accordance with standard dental medicine practices, dentist must clean medical device before inserting it into the oral cavity. The dentist must provide the patient with instructions for adequate cleaning of the medical device and follow up the patient at regular appointments to assess the condition of the medical device. We recommend cleaning with soft brush and non-abrasive toothpaste. Non-compliance with these instructions may result in increased biofilm formation or discoloration of the denture base material.

General warnings for the dental technician: Do not use material after expiry date. Expiry date is printed on the packaging. Keep out of reach of children. Keep the liquid and powder in the original packaging only. Do not use if the primary packaging is damaged. Avoid direct contact between skin and uncured material as well as inhalation of the vapours. It is necessary to wash your hands after handling the uncured material. Close containers immediately after use. Ensure efficient ventilation of the workspace. Ensure sufficient extraction of the particles from the workplace during the finishing steps (trimming, grinding, polishing). Wear suitable personal protective equipment (PPE) during all procedures. If no technical protective measures are taken, irritation of the eyes, skin and respiratory organs may occur during the handling with uncured material and finishing procedures. Consult the Safety Data Sheet (SDS) for proper PPE. (available on the website https://www.polident.si). Do not store the powder at temperatures above 50°C for long period of time.

Warnings for liquid component: Liquid is highly flammable. Irritating to eyes, respiratory system and skin. May cause sensitization by skin contact (allergic skin reaction). Keep away from heat, sparks, open flames, and other sources of ignition. Do not smoke. Avoid breathing vapours. Bottles must be tightly closed and stored in well-ventilated place. Do not get it on eyes, on skin, or on clothing. Wear protective gloves/protective clothing/eye protection. Take off immediately all contaminated clothing. Rinse skin with water or shower. Do not eat/drink. Take precautionary measures against static discharges.

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