ORGANIC PMMA FD

Instructions for use

EN

INTENDED PURPOSE:

ORGANIC PMMA FD disc is suitable for milling full and partial dentures and temporary immediate load implant supported dentures.

COMPOSITION: PMMA material with embedded pigments and residual methyl methacrylate less than 1%

INDICATION: Edentulism, partial edentulism

CONTRAINDICATIONS: Do not use the product in case of known allergies or hypersensitivity to acrylic materials.

WEARING TIME: The recommended wearing time is up to 5 years.

INSTRUCTIONS FOR USE

NOTE! Before use of the disc carefully read the label of the disc to choose the correct height of the disc and shade.

Design rules to consider within CAD software for the design of denture bases

Design the denture base with suitable CAD dental software.

Fix the disc in the milling machine holder. Pay attention to the orientation of the multilayer disc. The lighter shade should be on top. Load the CAM -file on the milling machine. Milling strategies for PMMA must be considered according to the milling machine instructions. The denture is separated from the disc with standard grinding tools for PMMA resin after milling.

Bonding of denture teeth

Bond the artificial teeth to the denture base with PMMA resin or light curing composite material according to manufacturer's instructions. The artificial teeth must be roughened and prepared according to manufacturer's instructions. The denture base does not require special treatment in case of PMMA bonding resin but usually requires a primer treatment in case of composite bonding material according to manufacturer's instructions. The denture can be finished with standard grinding and polishing tools for PMMA resin after the bonding of artificial teeth.

Finishing

Tools for acrylate are used for finishing and final adjustments. For polishing use grinding paper, polishing burrs, pumice stones, pastes and soft polishing brushes suitable for acrylic resins to achieve glossy surface.

CLEANING: Make sure to thoroughly clean the restoration to remove any residue of the milling and finishing processes using common dental laboratory practice (e.g. steam jet or mild soap/detergent).

INFORMATION FOR DENTIST: In accordance with standard dental medicine practices, dentist must clean medical device before inserting it into the oral cavity. The dentist must provide the patient with instructions for adequate cleaning of the medical device and follow up the patient at regular appointments to assess the condition of the medical device. We recommend cleaning with soft brushes and mild soap for removable medical device and non-abrasive tooth pastes for fixed restorations.

STORAGE: Keep away from direct sunlight!

DISPOSAL: Disposal should be made in accordance with local, state and national legislation.

WARNINGS: Professional use only! Dental CAD/CAM skills are required to properly process the product. The discs must be milled using a compatible CAD/CAM system.

It is important not to overheat the material!

Do not use the material after expiry date. Expiry date is printed on the packaging.

PMMA dust may develop during the milling process, leading to mechanical irritation of the eyes, skin, and respiratory tract. Make sure to use the extraction system in the workplace. When processing PMMA discs, wear personal protective equipment (protection mask for dust, safety glasses, ...).

To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

Further safety-relevant information can be found in the material safety data sheet available on request.

In case that any serious incident has occurred in relation to the device it should be reported to the manufacturer (Polident d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, where the user is established.

Technical Parameters

	Specification	Typical average value	Method according to
Flexural Strength	min. 65 MPa	90 MPa	ISO 20795-1:2013
Flexural modulus	min. 2000 MPa	2700 MPa	ISO 20795-1:2013
Water sorption	< 32 µg/mm ³	25 µg/mm ³	ISO 20795-1:2013
Water solubility	< 1,6 µg/mm ³	0,6 µg/mm ³	ISO 20795-1:2013
Residual methyl methacrylate	< 2,2%	0,5 %	ISO 20795-1:2013

List of symbols

catalogue number

keep away from sunlight

batch code

Medical Device

use by date

CE marking

consult instructions for use

Unique Device Identification

manufacturer

distributor

Organical CAD/CAM GmbH
Ruwersteig 43, 12681 Berlin, Germany
<http://www.organical-cadcam.com>

POLIDENT d.o.o.
Volčja Draga 42, 5293 Volčja Draga, Slovenija
<http://www.polident.si>

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