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DD poly Y ML

	Spezifikation	Typischer Durch-schnittswert	Methode gemäß
Biegebruchfestigkeit	> 65 MPa // > 50 MPa	90 MPa // 120 MPa	EN ISO 20795-2 // EN ISO 10477
Biegemodul	> 2000 MPa	2700 MPa	EN ISO 20795-2
Sorption	< 32 µg/mm³	25 µg/mm³	EN ISO 20795-2
Löslichkeit	< 1,6 µg/mm²	0,6 µg/mm²	EN ISO 20795-2
Restmonomer	< 2,2%	0,5 %	EN ISO 20795-2
Vickers-Härte	> 24 HV	26 HV	/

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Dimensions: Shape: 98,5 mm - shoulder Thickness: 16, 20 mm

Shades: OM2 - OM3 - A1 - A2 - A3 - A3.5 - A4 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4

INTENDED PURPOSE:

DD poly Y ML disc is suitable for milling:

- long-term temporary crowns and bridges,
- individual denture teeth, full arch and tooth segments to manufacture denture.

Multilayer **DD poly Y ML** disc is composed by several layers of transient color. The layering concept is optimized to match the natural color gradient of natural teeth. The final result is a highly aesthetic and economical restoration.

COMPOSITION: PMMA material with embedded pigments and residual methyl methacrylate less than 1%

WEARING TIME: The recommended wearing time is up to 5 years

CONTRAINDICATIONS: Do not use the product in case of known allergies or hypersensitivity to acrylic materials.

NOTE! Before use of the disc carefully read the label of the disc to choose the correct height of the disc and shade.

INTENDED USE: CROWNS AND BRIDGES

Indication: damaged teeth; partial edentulism

NOTE! Cantilever design of bridges and application on implants are not recommended because of the increased risk of fracture. If implant restoration is required, the screw channel may further increase the risk of fractures.

Rules to consider within the design of C & B:

- Keep a minimum wall thickness of 1,0 mm in occlusal region and 0,6 mm in cervical region.
- The connector section should be more than 10 mm² in anterior region and more than 15 mm² in posterior region.
- A max of two pontics is allowed between two abutment teeth in posterior region.
- Pay attention to the positioning of the restoration within the MULTILAYER disc. The enamel layer is marked on the top of the disc.

For finishing processes see FINISHING STEP

Cementing in mouth

Use adequate cement for cementing the crown or bridge. Bonding surface of crowns and bridges should be cleaned and dried before cementing. Follow the instructions provided by the manufacturer of dental cement.

INTENDED USE: DENTURE TEETH

Indication: Edentulism, partial edentulism

NOTE! Application on implants is not recommended due to the increased risk of fractures. When implant restoration is required, the screw channel may further increase the risk of fractures.

Rules to consider within the design of denture teeth:

- Design the teeth, set of teeth or full arches of teeth within the denture design software.

- An optimal design thickness in the vertical dimension should be maintained.

- Mill the teeth on a CNC milling machine in accordance with the machine producer instructions.

In case of multilayer discs the incisal side should be on top corresponding to the incisal or occlusal direction of teeth.

- Use suitable tools for acrylate to separate the teeth from the disc, and for finishing and final adjustments.

Bonding to PMMA denture base:

Roughen and wet the teeth and base with monomer. Apply adequate bonding agent in the sockets and press the corresponding teeth in position. The bonding agent must cover the entire tooth surface when pressed on the denture base. Clean the residual resin.

FINISHING STEP:

- Tools for acrylate are used for finishing and final adjustments. For polishing use pumice stone, pastes and soft polishing brushes to achieve a glossy surface.

CLEANING: Make sure to thoroughly clean the restoration to remove any residue of the milling and finishing processes using common dental laboratory practice (e.g. steam jet or mild soap/detergent).

INFORMATION FOR DENTIST: In accordance with standard dental medicine practices, dentist must clean medical device before inserting it into the oral cavity. The dentist must provide the patient with instructions for adequate cleaning of the medical device and follow up the patient at regular appointments to assess the condition of the medical device. We recommend cleaning with soft brushes and mild soap for removable medical device and non-abrasive tooth pastes for fixed restorations.

STORAGE: Keep away from direct sunlight!

DISPOSAL: Disposal should be made in accordance with local, state and national legislation.

CAUTION! For US users only: Federal law restricts this device to sale by or on the order of a dentist.

WARNINGS: Professional use only! Dental CAD/CAM skills are required to properly process the product. The discs must be milled using a compatible CAD/CAM system.

It is important not to overheat the material!

Do not use the material after expiry date. Expiry date is printed on the packaging.

PMMA dust may develop during the milling process, leading to mechanical irritation of the eyes, skin, and respiratory tract. Make sure to use the extraction system in the workplace. When processing PMMA discs, wear personal protective equipment (protection mask for dust, safety glasses, ...).

To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

Further safety-relevant information can be found in the material safety data sheet available on request.

In case that any serious incident has occurred in relation to the device it should be reported to the manufacturer (POLIDENT d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, where the user is established.

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Dental Direkt GmbH Industriezentrum 106–108 32139 Spenge | Germany **www.dentaldirekt.de**